



Cellometer[®] Vision CBA

Image Cytometry System for Cell-Based Assays

Specification Sheet

Includes:

- Cellometer Vision HSL Instrument
- Cellometer Software
- FCS Express 4 analysis software (2 licenses)
- Two Fluorescent Optics Modules
- Power Supply
- USB 2.0 Connection Cable
- Phone/online applications support during set-up

Magnification Selection:

- 5x (ex: most mammalian cells)
- 10x (ex: yeast, platelets)

Imaging

Performance:

	5x magnification	10x magnification
Cell Size	5 – 300* microns	0.5 – 25 microns
Conc. Range:	10 ⁵ – 10 ⁷ cells/ml	10 ⁵ – 10 ⁷ cells/ml

Bright field imaging, fluorescent imaging and pattern-recognition software to quickly and accurately decluster, identify and count individual cells.

* Cellometer CHT4-PD300 Slides are required for cells > 80 microns in diameter

Instrument Specifications:

Weight:	25.0 lbs. (11.3 kg)
Dimensions:	Width: 6.0" (15.2 cm) Depth: 8.5" (21.6 cm) Height: 14.0" (35.6 cm)
Input to Power Adapter:	100-240 VAC, 50/60 Hz, 1.0A
Output to Instrument:	12 VDC, 3.34A

PC / Laptop

Requirements:

(If purchasing

Cellometer

without PC laptop)

- Windows Vista, Windows 7 or Windows 8
- 1.6 GHz or higher processor
- 2GB RAM
- USB 2.0 port
- Display resolution: 1024 x 768 pixels or higher

Intel Mac with Parallels/Windows



Nexcelom Bioscience LLC. | 360 Merrimack Street, Building 9 | Lawrence, MA 01843
Telephone: 978.327.5340 | Fax: 978.327.5341 | Email: info@nexcelom.com | www.nexcelom.com
ISO 9001:2008 registered

1001092 Rev. D



Cellometer[®] Vision CBA

Image Cytometry System for Cell-Based Assays

Specification Sheet

**Available
Fluorescent Optics
Modules:**

VB-450-302

Excitation / Emission: 375nm/450nm

Example Fluorophores:

- DAPI
- Hoechst 33342 & 33258

VB-535-402

Excitation / Emission: 470nm/535nm

Example Fluorophores:

- Acridine Orange (+DNA)
- GFP
- FITC

VB-595-502

Excitation / Emission: 525nm/595nm

Example Fluorophores:

- Propidium Iodide (PI)
- Ethidium Bromide (EB)
- R-phycoerythrin (PE)

VB-660-502

Excitation / Emission: 540nm/660nm

Example Fluorophores:

- Propidium Iodide (PI)
- Ethidium Bromide (EB)
- 7-AAD

VB-695-602

Excitation / Emission: 630nm/695nm

Example Fluorophores:

- Allophycocyanin (APC)
- Cy5[®]

For more information, email support@nexcelom.com or call: 978-327-5340.



Nexcelom Bioscience LLC. | 360 Merrimack Street, Building 9 | Lawrence, MA 01843
Telephone: 978.327.5340 | Fax: 978.327.5341 | Email: info@nexcelom.com | www.nexcelom.com
ISO 9001:2008 registered

1001092 Rev. D