

Cellometer® Mini

Cell Counter

➔ Affordable Automated Cell Counting

➔ “Wicked Fast” Results

Concentration and Viability in < 10 seconds!


➔ Improved Accuracy vs. Manual Counting

Superior identification of live vs. dead cells, individual cells in clusters, and irregular-shaped cells


➔ Improved Precision vs. Manual Counting

- Eliminate manual counting and recording errors
- Reduce inter-operator variability
- Eliminate wash steps and potential cross-contamination

➔ View, Print, and Save Counted Images, Data Reports, and Cell Size Histograms



Results	
Assay: Cell line, viability trypan blue	
Sample ID: Cell Line Dilution Factor: 2.00	
Count	Concentration
Total: 155 cells	2.20x10 ⁶ cells/mL
Live: 148 cells	2.10x10 ⁶ cells/mL
Dead: 7 cells	9.92x10 ⁴ cells/mL
Viability: 95.5%	



For more information on the new Cellometer Mini, contact Nexcelom at 978-327-5340 or info@nexcelom.com!



w: www.nexcelom.com
p: 978.327.5340
e: info@nexcelom.com

Cellometer®
Simply Counted

Celigo®
Image Cytometer

Nexcelom products are for RESEARCH USE ONLY and are not approved for diagnostic or therapeutic use.
© Copyright 2014 Nexcelom Bioscience LLC. All Rights Reserved.

Nexcelom
Bioscience

1001201 Rev.B 07/14

